

Luminescence Assay for Alkaline Phosphatase Activity of C3H10T(1/2) Cells

Reagents

Luminometer assay plate

Lysis buffer (100 mM Tris pH 9.5, 250 mM NaCl, 25 mM MgCl₂, 1% Triton X-100)

CDP-Star Chemiluminescence reagent (PerkinElmer)

Procedure

- 1) Seed 96-well tissue culture plate with 5000 C3H10T1/2 cells/well in DMEM containing 10% fetal bovine serum and penicillin/streptomycin.
- 2) Allow cells to grow to confluency (typically 36-48 hours).
- 3) Switch cells to DMEM + 0.5% calf serum medium containing appropriate signaling agonists and or antagonists and controls. Continue growth of cells for 36-48 hrs.
- 4) Remove media and add 50 µl of lysis buffer to each well. Incubate with rocking for 45 minutes.
- 5) Transfer 10 µL of lysate to 96-well assay plate.
- 6) Add 50 µL of room temperature CDP-Star reagent to each lysate-containing well of the assay plate, and incubate for 15 minutes at room temperature in the dark.
- 7) Read plate on luminometer. One second or less integration time should be sufficient.

Notes

- 1) Cells may be seeded at higher density to reduce time to confluency.
- 2) Amount of lysate used in assay, as well as, amount of reagent are flexible but should not exceed a 1:2 ratio.
- 3) Incubation time before plate reading may have subtle effects on data and should be held constant for best comparison between experiments.