

cDNA Synthesis and Amplification for Microarrays

This procedure was adapted from protocols for TRIzol (Invitrogen) and the WTA2 TransPlex Whole Transcriptome Amplification kit (Sigma-Aldrich).

Materials

Sigma WTA2 kit, 100 reactions
QIAquick PCR cleanup kit
Agilent DNA 1000 kit for bioanalyzer
Linear acrylamide (5 mg/mL), Ambion
RNase-free chloroform
RNase-free 75% ethanol
RNase-free isopropanol
RNase-free microcentrifuge tubes and pipette tips
Siliconized 1.5 mL microcentrifuge tubes

RNA isolation

- 1) To 1 mL of thawed RNA-containing TRIzol solution, add 5 μ L linear acrylamide and mix tube.
- 2) Add 200 μ L of chloroform and shake vigorously for 20 seconds. Let sit 2-3 min. room temperature
- 3) Centrifuge at 12,000 x g for 15 min in at 4 °C.
- 4) Transfer aqueous phase to a new siliconized 1.5 mL microcentrifuge tube with a P200 and add 0.5 mL of isopropanol. Mix by inverting the tube a few times, and incubate the solution for 10 min at room temperature. Invert the tube several times and centrifuge at 12,000 x g for 15 min at 4 °C.
- 5) Using a P1000, carefully remove supernatant. Watch for the pellet, which will be slightly yellow in color at this step.
- 6) Add 1 mL of 75% ethanol to the pellet, which should make it appear white. Invert the tube several times, and then centrifuge at 7,500 x g for 5 min at 4 °C.
- 7) Carefully remove as much supernatant as possible with a P1000. Briefly centrifuge the tube and remove the remaining supernatant with a P10. Air dry briefly but do not allow the pellet to dry completely. Resuspend the pellet in 2.8 μ L of RNase-free H₂O.
- 8) Transfer this solution to an RNase-free PCR tube.

Reverse Transcriptase reaction

- 1) Thaw the following WTA2 kit components:
 - a) Library synthesis buffer
 - b) Library synthesis solution
 - c) Library synthesis enzyme
 - d) Nuclease-free water.

Mix library synthesis buffer and library synthesis solutions and dissolve any precipitates by briefly heating to 37 °C with thorough mixing.

- 2) Add 0.5 μL of library synthesis solution to the 2.8 μL solution of RNA.
- 3) Prime the RNA in a thermocycler: 70 $^{\circ}\text{C}$ for 5 min then hold 18 $^{\circ}\text{C}$. While priming, make library synthesis mix (volumes/reaction):
 - a) 0.5 μL library synthesis buffer.
 - b) 0.78 μL water.
 - c) 0.5 μL library synthesis enzyme.
- 4) Add 1.78 μL of the library synthesis mix to the primed RNA at 18 $^{\circ}\text{C}$. There should be 5.08 μL of total solution per tube.
- 5) Incubate this mixture in the thermocycler:
 - a) 18 $^{\circ}\text{C}$ for 10 min
 - b) 25 $^{\circ}\text{C}$ for 10 min
 - c) 37 $^{\circ}\text{C}$ for 30 min
 - d) 42 $^{\circ}\text{C}$ for 10 min
 - e) 70 $^{\circ}\text{C}$ for 20 min
 - f) Hold 4 $^{\circ}\text{C}$

Amplification reaction

- 1) Thaw the amplification mix and 10 mM dNTP mix.
 - a) Prepare the following master mix (given volumes are for one amplification reaction): 63.5 μL RNase free water
 - b) 8 μL amplification mix
 - c) 1.6 μL 10 mM dNTP mix
 - d) 0.8 μL 1:1000 SYBR green (optional for qPCR, if omitted use equal volume of water)
 - e) 0.8 μL amplification enzyme
- 3) Add 70 μL of master mix to each cDNA sample (will have 75.1 μL total volume)
- 4) Incubate in thermal cycler using the following parameters:
 - a) 94 $^{\circ}\text{C}$ for 2 min
 - b) 22 cycles of 94 $^{\circ}\text{C}$ for 30 sec, 70 $^{\circ}\text{C}$ for 5 min
 - c) If monitoring by qPCR, proceed 2-3 cycles after plateau.

Purify using QIAquick PCR cleanup kit

- 1) Elute in 30 μL EB according to the manufacturer's instructions.
- 2) Quantify DNA concentrations on a Nanodrop spectrophotometer
Should be $>200 \text{ ng}/\mu\text{L}$ ($> 10 \mu\text{g}$)
- 3) Run sample on Agilent Bioanalyzer using DNA 1000 kit (1.0 μL)
Should have an average cDNA length of approximately 400 bp