

High-Sensitivity Silver Staining of SDS-PAGE Gels

Reagents

- 50% MeOH, 12% acetic acid, and 100 μ L 37% paraformaldehyde solution per 200 mL (Fix Solution)
- 50% MeOH (Wash Solution)
- 0.02% (w/v) $\text{Na}_2\text{S}_2\text{O}_3$ (Pretreatment Solution)
- 0.2% (w/v) AgNO_3 (Silver Solution)
- 6% Na_2CO_3 and 4 mL 0.02% $\text{Na}_2\text{S}_2\text{O}_3$ /100 μ L 37% paraformaldehyde solution per 200 mL (Develop Solution)
- 50% MeOH, 12 % acetic acid (Stop Solution)
- 30% MeOH, 5% glycerol (Drying Solution)

Procedure

- 1) Treat gel with the Fix Solution for at least 1 hour (can be done overnight). 200 mL of Fix Solution is enough for 2-4 minigels.
- 2) Wash the gel with the Wash Solution for 3 x 8 minutes.
- 3) Treat the gel with the Pretreatment Solution for *exactly* 1 minute. The gel will float at first, and so care should be taken to uniformly treat the gel by agitation.
- 4) Rinse the gel with MilliQ water for *exactly* 3 x 20 seconds.
- 5) Incubate the gel with the Silver Solution for 20 minutes
- 6) Rinse the gel with MilliQ water for *exactly* 2 x 30 seconds.
- 7) Incubate the gel with the Develop Solution with agitation until the desired level of staining is achieved (typically 2-10 minutes). Be sure to have the Stop Solution ready.
- 8) Stop the staining by treating the gel with the Stop Solution for at least 10 minutes.
- 9) If you want to preserve the gel between sheets of porous plastic, wash the gel with the Wash Solution for 2 x 10 minutes. Then incubate the gel with the Drying Solution for a maximum of 1 hour (the glycerol will begin to leach out the silver after 1-2 hours).

Notes

- 1) This procedure can reliably detect protein levels down to about 5 ng.
- 2) All steps should be conducted with mechanical or manual swirling.
- 3) For more control over the staining rate, the Develop Solution can be removed from the gel when about 80% of the desired intensity is achieved. The gel will continue to develop slowly with the remaining reagents trapped in the gel.