

## 2 $\theta$ – $\omega$ scan using PDS and PIXcel in Bragg-Brentano geometry (X'Pert\_2 only)

### I. Login

1. Start **Data Collector**.
2. Type your *UserName* and *Password*.
3. Select *Instrument – Connect*.
4. Choose Configuration **PDS + PIXcel**.
5. Click *OK*.

### II. Hardware Setup

1. X-ray Tube is in *Line Focus*.
2. Goniometer Resolution set to “Normal 0.001 deg”.
3. Incident Beam Optics – **PDS (Programmable Divergence Slit)**  
*Note: if you have to change incident beam optics please first turn Automatic attenuator to “Activate” status and then unplug attenuator cable.*
  - a. If needed, insert Ni filter.
  - b. If sample is smaller than 25 mm, insert correct size *Mask*.
4. Diffracted Beam Optics – **PIXcel**
  - a. Insert Soller Slit.

### III. Data Collector Software

1. Select the *Incident Beam* optics tab
  - a. Double click any item. Incident beam optics window will appear.
  - b. Go through all tabs and select proper optic components.
  - c. *Divergence Slit* – “Fixed” mode and slit set to  $1/32^\circ$ .
  - d. Set *Automatic Attenuator Usage* to “Do not switch” and “Activated”.
2. Select the *Diffracted Beam* optics tab
  - a. Double click any item. Diffracted beam optics window will appear.
  - b. Go through all tabs and select proper optic components.
  - c. *Detector* – “Receiving Slit” mode and active length set to 0.16 mm.
3. Select *Instrument Settings* tab.
  - a. Double click any item in the tree view to prompt another window.
  - b. Press *X-ray* tab. Set generator power to 45 kV and 40 mA.

### IV. Diffractometer Zero Alignment

1. In *Instruments Settings* check Z position. If it is larger than 5 mm move it back to at least 5 mm.
2. Move all other motors to zero positions.
3. From Menu Select *Measure – Manual Scan*.
4. From the *Scan Axis* drop down menu select 2Theta.
5. Enter *Range*  $1^\circ$ , *Step Size*  $0.005^\circ$ , and *Time per Step* 0.1sec. Then press *Start*.
6. After scan is finished, right click on mouse and select *Move Mode*. Move Scan Axis to the center of the mass of the peak.
7. In *User Settings – Sample Offsets* set current 2Theta position to zero.

### V. Sample Mounting

1. Mount the sample.
2. If in the *Instrument Settings* tab X = 0.0 and Y = 0.0, beam is positioned at the center of a sample stage (aluminum disk).

### VI. Moving Sample into the Beam Position

1. Using supplied Micrometer

- a. Mount micrometer onto the MRD cradle. Close the doors.
  - b. In the *Instruments Settings* window, double click any item in the tree view to prompt another window.
  - c. Move Z until micrometer reads  $(1.00 \pm 0.02)$  mm. This is the correct sample height.
2. Using direct beam
    - a. Note the direct beam intensity. In *Instruments Settings* move Z to higher values until intensity starts to drop.
    - b. In *Manual Scan* window from the *Scan Axis* drop down menu select Z. Enter *Range* 2 mm, *Step Size* 0.01 mm, and *Time per Step* 0.2 sec. Then press Start.
    - c. After scan is finished, right click on mouse and select *Move Mode*. Move *Scan Axis* to the intensity value corresponding to  $\frac{1}{2}$  of the direct beam intensity.
    - d. Z is aligned. Close shutter.
 

*Note: sample height alignment using direct beam will work correctly if sample size is larger than the beam size in axial direction.*

## VII. Measurement

1. Select *Incident Beam* optics tab.
2. The *Divergence Slit* can be used in “Fixed” and “Automatic” modes. If “Fixed” mode is set, select the divergence of the slit to desired value from  $1/32^\circ$  to  $4^\circ$ . If “Automatic” mode is set, select a proper irradiated length of the sample.
3. Set *Automatic Attenuator Usage* to “Do not switch” and “Deactivated”.
4. Click on *Diffacted Beam* optics tab.
5. Select *Detector* tab. Set detector mode to “Scanning” and active detector angle set to  $2.511^\circ$ .
6. Simplest way to execute scan is to do a *Manual Scan*. It is a relative scan i.e. executed around current goniometer position with the range specified in *Manual Scan* window.
7. To do *2Theta – Omega* scan first move 2Theta and Omega to middle positions of the scan range.
8. In *Manual Scan* window select *Scan Axis* “2Theta – Omega” and appropriate *Range*, *Step Size* and *Time per Step*.
9. Click *Start*. When scan is completed, save it through *File – Save As* menu. Manual Scan will be lost if it is not saved.

*Note: In symmetrical scan mode Omega axis is always  $\frac{1}{2}$  of the 2Theta axis.*

## VIII. Logging out

1. Close the shutter.
2. Switch *Automatic Attenuator* to “Do Not Switch” and “Activated”.
3. Move all angles to zero positions and Z to 5 mm.
4. Lower the power of the x-ray tube to 40 kV and 10 mA.
5. Close *Data Collector*.