

## Miscellaneous Solutions

### 1x WB 1x Ni-NTA Chromatography Loading/Wash Buffer, pH 8.0

Chemical	FW	to add for 1 l
300 mM NaCl	58.44	17.5 g
20 mM HEPES	238.3	4.77 g
20 mM Imidazole	68.08	1.36 g

Bring to 950 ml with H<sub>2</sub>O, bring to pH 8.0 (~35 drops 10 N NaOH), then bring to 1 l.

### 10x WB 10x Ni-NTA Chromatography Loading/Wash Buffer, pH 8.0

Chemical	FW	to add for 1 l
3 M NaCl	58.44	175.3 g
200 mM HEPES	238.3	47.7 g
200 mM Imidazole	68.08	13.6 g

Bring to 950 ml with H<sub>2</sub>O, bring to pH 8.0 (~100 drops 10 N NaOH), then bring to 1 l.

### 1x EB Ni-NTA Chromatography Elution Buffer, pH 8.0

Chemical	FW	to add for 1 l
300 mM NaCl	58.44	17.5 g
20 mM HEPES	238.3	4.77 g
250 mM Imidazole	68.08	17.0 g

Bring to 975 ml with H<sub>2</sub>O, bring to pH 8.0 (~50 drops 37% HCl), then bring to 1 l.

### Crystallization Buffer 10 mM NaCl, 10 mM HEPES, pH 7.4

Chemical	FW	to add for 1 l
10 mM NaCl	58.44	0.58 g
10 mM HEPES	238.3	2.38 g

Bring to 975 ml with H<sub>2</sub>O, bring to pH 7.4 (~7 drops 10 N NaOH), then bring to 1 l.

**NMR Buffer** 50 mM NaCl, 10 mM Phosphate pH 7.4, 0.02% NaN<sub>3</sub>

Chemical	FW	to add for 1 l
50 mM NaCl	58.44	2.92 g
2.3 mM Monosodium phosphate, monohydrate	137.99	0.31 g
7.7 mM Disodium phosphate, heptahydrate	268.07	2.07 g
0.02% NaN <sub>3</sub> (Sodium Azide)	65.01	0.2 g

Bring to 1 l.

**Coomassie Destain #1** Rinse Solution

Chemical	to add for 1 l
95% Ethanol	280 ml
Glacial Acetic Acid	100 ml
Water	620 ml

**Coomassie Destain #2** Drying Solution

Chemical	to add for 1 l
Glycerol	100 ml
Glacial Acetic Acid	50 ml
Water	850 ml

**Anion Exchange Buffer A** 10 mM NaCl, 30 mM Tris, pH 8.0

Chemical	FW	to add for 1 l	
		4°C	25°C
10 mM NaCl	58.44	0.58 g	0.58 g
30 mM Tris			
Tris-HCl	157.6	3.97 g	2.67 g
Tris (Base)	121.1	0.58 g	1.59 g

**Anion Exchange Buffer B** 1 M NaCl, 30 mM Tris, pH 8.0

Chemical	FW	to add for 1 l	
		4°C	25°C
1 M NaCl	58.44	58.4 g	58.4 g
30 mM Tris			
Tris-HCl	157.6	3.97 g	2.67 g
Tris (Base)	121.1	0.58 g	1.59 g

**1M Tris, pH 8** Will be pH 8 when diluted to 30 mM at the given temperature

Chemical	FW	to add for 1 l	
		4°C	25°C
1 M Tris			
Tris-HCl	157.6	132.2 g	88.8 g
Tris (Base)	121.1	19.4 g	53.0 g

**5x M9 Salts** - Based on Maniatis, A2.2

Chemical	FW	to add for 1 l
Na <sub>2</sub> HPO <sub>4</sub> · 7H <sub>2</sub> O	268.07	64 g
KH <sub>2</sub> PO <sub>4</sub>	136.09	15 g
NaCl	58.44	2.5 g
NH <sub>4</sub> Cl (leave out for isotopic labeling)	53.50	5.0 g

Sterilize by autoclaving.

**1 M MgSO<sub>4</sub>** - Based on Maniatis, A2.2

Chemical	FW	to add for 100 ml
MgSO <sub>4</sub>	120.37	12.0 g

Sterilize by autoclaving.

**1 M CaCl<sub>2</sub>** - Based on Maniatis, A2.2

Chemical	FW	to add for 100 ml
CaCl <sub>2</sub> · H <sub>2</sub> O	147.02	14.7 g

Sterilize by autoclaving.

**M9 Minimal Media** - Based on Maniatis, A2.2

Ingredient	to add for 1 l
H <sub>2</sub> O (Sterile)	775 ml
5x M9 Salts (Sterile)	200 ml
1 M MgSO <sub>4</sub> (Sterile)	2 ml
20% Glucose Solution (Filter Sterilized)	20 ml
1 M CaCl <sub>2</sub> (Sterile)	0.1 ml
1000x Antibiotic (Filter Sterilized)	1 ml
if cells are $\Delta$ leu	
20 mg/ml L-Leucine (Filter Sterilized)	2 ml

To make M9 agar plates, make agar at 20 g/l and autoclave. Then make media as above, but use the 20 g/l agar instead of H<sub>2</sub>O to obtain a final agar concentration of 15 g/l.

**1000X DNase/RNase stock :**

Ribonuclease A (Bovine Pancreas, Sigma R-5500)	10 mg
Deoxyribonuclease I (Bovine Pancreas, Sigma D-4527, 2000 units/mg)	5 mg
PBS	1 ml

0.2  $\mu$ m filter sterilize.